

# Large complexes: Cloning strategy, production and purification

Eric Durand, Roland Lloubès

### ► To cite this version:

Eric Durand, Roland Lloubès. Large complexes: Cloning strategy, production and purification. Methods in Molecular Biology: Bacterial protein secretion systems, 2017, 10.1007/978-1-4939-7033-9\_24. hal-01788569

## HAL Id: hal-01788569 https://amu.hal.science/hal-01788569

Submitted on 9 May 2018  $\,$ 

**HAL** is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers. L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.

#### 1Chapter 24. Large complexes: Cloning strategy, production and purification

2Running Head: Membrane protein complex purification

3Eric Durand and Roland Lloubes

4attention abbreviations

5Laboratoire d'Ingénierie des Systèmes Macromoléculaires UMR7255, Institut de Microbiologie 6de la Méditerranée, CNRS/Aix-Marseille Université, 31 chemin Joseph Aiguier, 13402 Marseille 7cedex 20, France 8 9E-mail: roland.lloubes@imm.cnrs.fr 10

#### 11i. Summary

12Membrane proteins can assemble and form complexes in the cell envelope. In Gram-negative 13bacteria, a number of multiprotein complexes including secretion systems, efflux pumps, 14molecular motors or pilus assembly machines comprise proteins from the inner and the outer 15membrane. Beside the structures of isolated soluble domains, only few atomic structures of these 16assembled molecular machines have been elucidated. To better understand the function and to 17solve the structure of protein complexes, it is thus necessary to design dedicated production and 18purification processes. Here, we present cloning procedures to overproduce membrane proteins 19into *E. coli* cells and describe the cloning and purification strategy for the Type VI secretion 20TssJLM membrane complex.

21

22ii. Key words: membrane protein complexes, *E. coli*, T7 overexpression, protein purification.23

#### 241. Introduction

25Protein overproduction results from cloning a gene of interest into a plasmid vector, downstream 26of a tightly regulated promoter, and from inducing its expression after plasmid transformation 27into a bacterial strain. For large protein complexes containing multiple subunits, the genes 28encoding the different subunits can be expressed under the control of an inducible promoter 29either from single plasmid containing a cluster of genes or from different compatible plasmids 30harboring single or multiple genes.

#### 311.1 Cloning vectors

32Several inducible promoters have been described and are available to overexpress a gene of 33 interest. These promoters are usually cloned into vectors that also contain the gene encoding the 34 cognate regulatory protein and a transcriptional terminator to prevent non-productive 35transcription from downstream gene (1-7). The tac and trc promoters that contain the -35 and -10 36sequences from the *trp* and *lacUV5* promoters, respectively, have been optimized for high 37 expression level (1). Other E. coli tightly regulated promoters such as the tetA promoter and its 38 repressor TetR, the *araBAD* promoter and its cognate activator AraC (3-6). In addition, 39heterologous combinations of promoter/regulator binding sequences can be used (2). The last 40 family of inducible promoters gathers sequences that are not recognized by the E. coli RNA 41polymerase but rather recognized by phage RNA polymerases such as the SP6, T3 and T7 42promoters. Vectors and strains that express the T7 RNA polymerase (T7RNAp) have been 43 extensively developed. Three independent methods are used to regulate the T7 expression 44systems. First, the expression of the chromosomally- or plasmid-encoded T7 RNAp gene can be 45itself under the control of an inducible promoter (see Note 1). Second, the inhibition of the T7 46RNAp basal activity can be controlled by producing the T7 lysozyme under constitutive or 47 regulated conditions (8, 9). Last, the transcription by the T7 RNAp can be repressed by the Lac 48 repressor, adding the *lac* operator sequence (10), the *lacI* gene being cloned on the expression 49vector.

51A large choice of cloning vectors with different characteristics is now available and can be 52selected based on the toxicity, stability and folding rate of the protein or protein complex to 53overproduce:

54- production of proteins in the cytoplasm, the periplasm or the membrane of *E. coli* (using the
55addition of synthetic N-terminal signal-sequences such as that of the OmpA and PelB proteins).
56- compatible T7 expression vectors containing one or two T7 promoters to overproduce protein
57complexes with up to eight subunits (see Duet Vectors from Novagen).

58- Protein tagging sequence (6 to 10 × His, *Strep*-Tag II...) or fusion partner (protein G, 59glutathione-S-transferase: GST, calmodulin-binding peptide: CBP, Maltose binding protein: 60MBP) to increase protein solubility or to simplify the purification using affinity chromatography 61techniques. In addition, these tag affinity sequences can be removed using specific protease 62(most used proteases: Tobacco Etch Virus Protease: TEV, enterokinase, Thrombin, Xa Factor, 63PreScission..). For this purpose, the corresponding protease recognition sequence is inserted 64either downstream (for N-terminal tagging) or upstream (for C-terminal tagging) of the tag 65sequence. After protease digestion, purified protein complex is submitted to a new purification 66step to remove the protease, the uncleaved protein and the tag peptide (for example affinity 67chromatography to remove the His-tag peptide and the His-tagged TEV).

68 Currently, cloning into these vectors is facilitated by PCR techniques based on gene and 69plasmid amplifications, leading to restriction site/ligation free cloning methods. These techniques 70have proven to be efficient to perform gene expression strategies in short delay (*see* section 3.1 71and **Note 2**).

If none of the overproduction strategies are sufficient to produce and purify protein
 73complexes, *in vitro* alternatives such as cell-free transcription-translation systems, based on

74T7RNAp and *E. coli* cell extracts, have been found efficient to produce milligram amounts of 75some membrane proteins *(11)*.

76

#### 771.2 Membrane protein complex overproduction and purification

78Membrane protein complexes have been successfully overproduced and purified based on pBAD 79or T7 vectors. For example, the Tol proteins from the Tol cell envelope molecular motor were 80overproduced and specifically radiolabeled with <sup>35</sup>S-Methionine after cloning into a pT7-driven 81vector using restriction sites (RS) and ligation techniques (RSI) (*12, 13, see* **Note 3 and Note 4**). 82Sub-complexes from the Ton system (*14-17*), the PomA-PomB flagellar rotor (*18*), the AcrAB-83TolC efflux pump (*19*), the Type IV and Type VI secretion systems (T4SS and T6SS, 84respectively; *20, 21*) and the β-barrel assembly machinery (*22*) have been also successfully 85purified using similar approaches.

In addition to the cloning strategy and the production to levels compatible with 87purification, studying these multiprotein machines also requires to extract and solubilize the 88protein complex without disrupting the contacts between the subunits. Finally, since the 89overproduction condition may induce stoichiometry artifacts, the tagging of the minor subunit or 90the specific tagging of different subunits should be performed to purify homogenous protein 91complexes (see below and section 3.2).

Here we describe the cloning strategy as well as the production, extraction, and 93purification of the T6SS TssJLM membrane complex (*21*).

94

952. Materials

#### 962.1 Cloning of the TssJLM complex

- 971. Desalted oligonucleotides of up to 35 bases, 5' ends unphosphorylated
- 982. Restriction enzymes (NdeI, XhoI, DpnI)
- 993. T4 DNA ligase
- 1004. Thermoblock (16°C-42°C)
- 1015. E. coli DH5a competent cells
- 1026. PCR thermocycler
- 1037. HPLC purified Mega-primer pairs ( $\geq$  50 bases), 5' ends unphosphorylated
- 1048. vector: pRSF-Duet1
- 1059. High fidelity Taq DNA polymerase (Pfu turbo, Agilent)
- 10610. dNTPs, 10 mM stock solution in water
- 10711. Agarose gel electrophoresis system

#### 1082.2 Production and purification of the TssJLM complex

- 1091. IPTG stock solution: 0.1 M in water
- 1102. BL21(DE3): E. coli strain B F<sup>-</sup> ompT gal dcm lon  $hsdS_B(r_B m_B) \lambda$ (DE3 [lacI lacUV5-T7 gene
- 1111 ind1 sam7 nin5]) [malB<sup>+</sup>] ( $\lambda$ <sup>S</sup>)
- 1123. Lysozyme stock solution: 10 mg/mL in water
- 1134. DNase stock solution: 10 mg/mL in water
- 1145. TCEP (Tris(2-carboxyethyl)phosphine hydrochloride) extemporaneously added
- 1156. EDTA stock solution: 0.5M in water
- 1167. MgCl<sub>2</sub> stock solution: 1 M in water
- 1178. Lysis buffer: 50 mM Tris-HCl, pH 8.0, 50 mM NaCl, 1 mM EDTA

1189. Solubilization buffer: 50 mM Tris-HCl, pH 8.0, 50 mM NaCl, 1 mM EDTA, 0.5% (w/v) n-

119dodecyl-β-D-maltopyranoside (DDM), 0.75% (w/v) decyl maltose neopentyl glycol (DM-NPG),

1200.5% (w/v) digitonin (Sigma-Aldrich)

12110. EDTA-free protease inhibitor tablets

122110. Affinity buffer: 50 mM Tris-HCl, pH 8.0, 50 mM NaCl, 0.05% (w/v) DM-NPG

12312. Imidazole-HCl, pH 8.0: stock solution 4M

12413. Incubation Shakers (16°-37°C)

12514. Centrifugation (from 5,000 to  $100,000 \times g$ )

12615. Glass potter

12716. Emulsiflex-C5 (Avestin)

12817. 5-mL StrepTrap HP and 5-ml HisTrap HP columns (GE Healthcare)

12918. Desthiobiotin

130

#### 1313. Methods

#### 1323.1: Cloning the TssJLM membrane complex

133To overproduce the T6SS TssJLM membrane complex, the *tssJ, tssL* and *tssM* genes have been 134assembled in an artificial operon with individual optimized ribosome-binding sites (*see* Figure 1351). In addition, based on previous data regarding permissive positions, each subunit has been 136tagged with a specific tag.

The pRSF–TssJ<sup>Strep</sup>-TssL<sup>FLAG</sup>-<sup>His6</sup>TssM plasmid was constructed by both RSl (*see* Note
and RSl-free (see Note 2, (23)) cloning methods (Figure 1).

1391. PCR amplify the *tssJ* gene (encoding the TssJ lipoprotein) using High fidelity Taq DNA140polymerase and primers containing the NdeI (5' primer) and the XhoI RS (3' primer) extensions.

141The *Strep*-Tag II DNA sequence followed by a stop codon was introduced in the 3' primer (Note1425). Check the correct amplification by agarose gel electrophoresis (AGE).

1432. Digest the PCR product and the pRSF-Duet vector with NdeI and XhoI restriction enzymes 144(sites present into the pRSF-Duet MCS2).3. Mix the digested PCR amplified fragment with 145about 50 ng of the digested pRSF-Duet vector, both digested with NdeI and XhoI, with a molar 146excess ratio of insert/vector of (2/1 to 5/1, inversely proportional to the DNA length of the 147insert).

1484. Add 1 U of T4 DNA ligase to the DNA mixture plus its specific buffer in a total volume of 15 149μl and incubate at 16-20°C for at least 2 hours.

1505. Transform *E. coli* DH5 $\alpha$  competent cells with about 30% of the ligation mix, plate on LB agar 151containing the appropriate antibiotic (in case of pRSF-Duet, use kanamycin 50 µg/mL) and 152incubate overnight at 37°C.

1536. Select positive clones containing the pRSF–TssJ<sup>Strep</sup> plasmid by colony PCR screening using 154the same primers. Extract the plasmid DNA and check the accuracy of the cloned sequence by 155DNA sequencing.

1567. PCR amplify the *tssL* gene using (i) a 5' primer containing a 5' extension (22bp) corresponding
157to the *Strep*-TagII and the stop codon followed by the RBS sequence (*see* Note 6) and the 5'158FLAG-*tssL* extension (63bp), and (ii) a 3' primer containing the 3'-*tssL* gene extension (35bp)
159including the stop codon followed by RBS, the ATG and the 5'-His6-*tssM* extension (14bp).
160Check correct amplification by AGE.

1618. PCR amplify the *tssM* gene using (i) the 5' primer complementary to the 3' end of *tssL* gene, 162RBS, His tag and 5' sequence of *tssM* gene (20bp) and (ii) the 3' primer containing the 3'-*tssM*  163gene extension (35bp) including the stop codon followed by the pRSF-Duet sequence extension164(35bp). Check correct amplification by AGE.

1659. Mix the *flag-tssL*, *his-tssM* PCR products and the pRSF–TssJ<sup>Strep</sup> plasmid and PCR amplify 166with High fidelity DNA polymerase in a single reaction as described for the RSI-free technique 167(*see* **Note 2**).

16810. Digest the mixture with *DpnI*. Transform into *E. coli* DH5 $\alpha$  competent cells and plate on LB 169agar containing the appropriate antibiotic and incubate overnight at 37°C.

17011. Select positive clones containing the pRSF–TssJ<sup>Strep\_FLAG</sup>TssL-<sup>6His</sup>TssM plasmid by colony 171PCR screening using the same primers. Extract the plasmid DNA and check the accuracy of the 172cloned sequence by DNA sequencing.

173

#### 1743.2: TssJLM membrane complex, extraction and purification

1751. Transform the expression vector (pRSF-TssJ<sup>Strep\_FLAG</sup>TssL-<sup>6His</sup>TssM) into the *E. coli* 176BL21(DE3) expression strain.

1772. Grow cells at 37 °C in 8 L of lysogeny broth (LB) to an optical density at 600 nm ( $A_{600}$ ) ~ 0.7. 178Induce the expression of the *tssJLM* genes with 1.0 mM IPTG for 16 h at 16 °C (*see* **Note 7**). 1793. Pellet cells by centrifugation at 7,000 × g for 20 min. Resuspend cell pellets in 300 mL of ice-180cold Lysis buffer supplemented with 1 mM TCEP, 100 µg/mL of DNase I, 100 µg/mL of 181lysozyme and with one tablet of EDTA-free protease inhibitor. Add MgCl<sub>2</sub> to the final 182concentration of 10 mM.

1834. Break the cell suspension with an Emulsiflex-C5 by 4 passages at 15,000 psi (100 MPa). 184Pellet unbroken cells by centrifugation at 7,000  $\times$  g (*see* **Note 8**).

1855. Pellet membranes by ultracentrifugation at  $98,000 \times g$  for 45 min.

1866. Resuspend membranes in 120 mL of Solubilization buffer supplemented with 1 mM TCEP at 18722°C and homogenize membranes mechanically with a potter (duration about 45 min) (*see* **Note** 188**9**).

1897. Clarify the membrane suspension by centrifugation at  $98,000 \times g$  for 20 min.

1908. Load the supernatant onto a 5-mL StrepTrap HP column and then wash with Affinity buffer at 1914 °C.

1929. Elute the TssJLM core complex in Affinity buffer supplemented with 2.5 mM desthiobiotin193into a 5-mL HisTrap HP column.

19410. Wash the HisTrap HP column in Affinity buffer supplemented with 20 mM imidazole and 195proceed to the elution of the TssJLM core complex in the same buffer supplemented with 500 196mM imidazole.

19711. Pool the peak fractions and load onto a Superose 6 10/300 column equilibrated in 50 mM 198Tris-HCl, pH 8.0, 50 mM NaCl, 0.025% (w/v) DM-NPG (*see* **Note 10**). The TssJLM complex 199elute as a single monodisperse peak close to the void volume of the column.

200

#### 2014. Notes

2021. *E. coli* T7RNAp expression systems result from (i) plasmid (pGP1-2) expression, the 203T7RNAp under the control of the Lambda P<sub>L</sub> promoter is regulated by the plasmid encoded 204temperature-sensitive C1-857 repressor (*24*); (ii) chromosomal expression: the T7RNAp gene is 205under the control of the *lacUV5* (*25*) or the *ara* promoter (*26*); (iii) besides tight control of the 206chromosomal T7RNAp expression upon AraC control (*26*), infection by phages (M13 mGP1-2 207or Lambda CE6) encoding the T7RNAp (*27, 28*) has proven to be efficient for the expression of 208toxic gene products.

2092. The-RSI-free cloning strategy allows the insertion of a gene at a precise position on a target 210plasmid but does not require the presence of RS. The restriction free cloning method consists in 211two sequential PCR amplifications using only two primers. These primers (with  $\geq 25$  bp 212homology extension) are designed to hybridize on the 5' and 3' ends of the gene of interest and 213contain additional extension to hybridize on the expression vector at a selected position 214downstream of the regulate promoter sequence. First, the gene of interest is amplified by PCR, 215then the amplified mega-primer pairs containing the gene are annealed to the vector of interest. A 216new amplification with High fidelity polymerase produce a linear gene-vector amplification (29). 217The gene is thus plasmid included into a nicked and circular DNA molecule. The PCR reaction is 218treated by DpnI RE to digest the unwanted Dam methylated plasmids templates. Competent 219E.coli cells are further transformed with the annealed DNA complex directly. An alternative-220method has been developed (30). It uses PCR amplified gene containing at least 15 bp extensions 221that are homolog to each end of the linearized vector (PCR amplified or RS digested). PCR DNA 222and vector treated with T4 polymerase, to create 5' overhangs, are annealed to form a 223recombinant plasmid.

224**3.** Addition of a T7 promoter by a RS-dependent cloning strategy. RS present in the upstream 225sequence of the gene of interest can be used to insert a synthetic T7 promoter DNA fragment 226which is formed by two overlapping oligonucleotides containing the T7 consensus promoter 227sequence of 23 bases corresponding to: 5'-taatacgactcactatagggaga-3'. The T7 promoter sequence 228is inserted upstream of the Ribosome Binding Site (RBS) of the gene of interest (it is also 229possible to insert the T7 promoter upstream of the natural promoter sequence). For this purpose, 230the synthetic DNA fragment contains additional 5' and 3' end extensions that are complementary 231of the sticky ends of the RS present in the plasmid (blunt ended RS are also possible to use but 232with lower ligation efficiency and random insertion). Two complementary desalted 233oligonucleotides (ON) of about 35 bases, 5' end unphosphorylated are hybridized following heat 234denaturation in ultrapure water and further cooling at room temperature. The annealed DNA 235fragment is further ligated into RS digested plasmid. It is possible to favor the selection of 236positive clones using Synthetic DNA that do not recreate the initial RS. Then, heat inactivated 237ligation mix is digested with the RE corresponding to the RS that has been destroyed upon DNA 238fragment ligation (Example: DNA fragment insertion into the EcoRI RS: after EcoRI digestion, 239the overhang EcoRI: 5'-AATTC... should be filled with the synthetic sequence: 5'-AATTX...., 240where the X nucleotide does not correspond to C nucleotide). It is noticeable that this fast 241technique does not require plasmid sequencing. It should be used for plasmids harboring 242resistance and regulatory genes in the opposite orientation from that of the T7 regulated gene(s) 243(*31*).

2444. The general RSI-dependent strategy is often used for cloning genes of interest into the MCS of 245expression vectors. The genes can be either purified from RS digestion *(32)* or obtained from 246PCR amplification using primers that contain additional RS extensions. Then, the PCR amplified 247DNA is RE digested and inserted into the cognate RS of the MCS present in the vector *(16)*. 2485. Since lipoproteins undergo a post-translational modification at their N-termini *(33)*, the *Strep*-249Tag II affinity tag sequence was introduced at the C-terminal of TssJ. It is important to note that 250the positions of affinity tags (Strep, FLAG and His Tags) have been rationally chosen to 251maintain functional proteins.

2526. To optimize the production of the TssL and TssM protein, the endogenous RBS-ATG 5' 253sequence in front of the genes were replaced (in the 5' primer sequence) by the consensual 254sequence "AAGGAGATATACAT*ATG*" (*34*) (RBS and start codon are presented, bold and 255italic letters respectively).

256 **7.** Growth and induction at 37°C for 3 hours led to very low biomass and protein yield. The 257culture was conducted at 37°C until  $A_{600} \sim 0.7$ -0.9 before shifting the incubator at 16°C prior to 258induction by IPTG. The final OD after the 16H induction was around 1.4-1.6. It is important not 259to induce around 0.4 OD since the cell would stop growing shortly after.

2608. It is important to note that whereas many cell disrupting protocol have been tried, only lysis 261with Emulsiflex produced a stable and homogenous sample of the TssJLM membrane core 262complex.

2639. This specific combination of detergents (same recipe than that used for the purification of the 264T4SS *(20)*) gave the higher extraction yield of the TssJLM membrane complex. Other isolated 265detergents (Triton-X-100, n-Octyl- $\beta$ -D-Glucoside, DDM or DM-NPG) gave poor extraction 266yields. All detergent buffers were prepared extemporally the day of the purification. Digitonin 267tends to precipitate in high salt buffers and without other detergents. Adding DDM and DM-268NPG prevent digitonin precipitation.

26910. DM-NPG was chosen because the solulilized TssJLM complex is stable and since this 270detergent is successful for high-resolution structural biology studies giving very clear and 271reproducible EM background after negative staining. Using low salt concentration (50 mM 272maximum) was key to prevent aggregation of the TssJLM complex. The affinity purification 273steps were of paramount importance to isolate a stoichiometric complex. Indeed, a far excess of 274the TssJ lipoprotein was eliminated during the second His column.

275

#### 2765. References

2771. Amann E, Brosius J, Ptashne M. (1983) Vectors bearing a hybrid trp-lac promoter useful for 278regulated expression of cloned genes in Escherichia coli. *Gene* 25, 167-178.

2792. Lutz R, Bujard H. (1997) Independent and tight regulation of transcriptional units in 280*Escherichia coli* via the LacR/O, the TetR/O and AraC/I1-I2 regulatory elements. *Nucleic Acids* 281*Res* 25,1203-1210.

2823. Cagnon C, Valverde V, Masson JM. (1991) A new family of sugar-inducible expression 283vectors for *Escherichia coli*. *Protein Eng* 4, 843-847.

2844. Skerra A. (1994) Use of the tetracycline promoter for the tightly regulated production of a 285murine antibody fragment in *Escherichia coli*. *Gene* 151,131-135.

2865. Guzman LM, Belin D, Carson MJ, Beckwith J. (1995) Tight regulation, modulation, and high-287level expression by vectors containing the arabinose PBAD promoter. *J Bacteriol* 177, 4121-2884130.

2896. Haldimann A, Daniels LL, Wanner BL. (1998) Use of new methods for construction of tightly 290regulated arabinose and rhamnose promoter fusions in studies of the *Escherichia coli* phosphate 291regulon. *J Bacteriol* 180, 1277-1286.

2927. Balzer S, Kucharova V, Megerle J, Lale R, Brautaset T, Valla S. (2013) A comparative 293analysis of the properties of regulated promoter systems commonly used for recombinant gene 294expression in *Escherichia coli*. *Microb Cell Fact*. 12, 26.

2958. Studier FW. (1991) Use of bacteriophage T7 lysozyme to improve an inducible T7 expression 296system. *J Mol Biol* 219, 37-44.

2979. Schlegel S1, Löfblom J, Lee C, Hjelm A, Klepsch M, Strous M, Drew D, Slotboom DJ, de 298Gier JW. (2012) Optimizing membrane protein overexpression in the *Escherichia coli* strain 299Lemo21(DE3). *J Mol Biol* 423, 648-659.

30010. Dubendorff JW, Studier FW. (1991) Controlling basal expression in an inducible T7301expression system by blocking the target T7 promoter with lac repressor. *J Mol Biol* 219, 45-59.

30211. Schwarz D, Junge F, Durst F, Frölich N, Schneider B, Reckel S, Sobhanifar S, Dötsch V,
303Bernhard F. (2007) Preparative scale expression of membrane proteins in *Escherichia coli*-based
304continuous exchange cell-free systems. *Nat Protoc* 2, 2945-2957.

30512. Guihard G, Boulanger P, Bénédetti H, Lloubes R, Besnard M, Letellier L. (1994) Colicin A 306and the Tol proteins involved in its translocation are preferentially located in the contact sites 307between the inner and outer membranes of *Escherichia coli* cells. *J Biol Chem* 269, 5874-5880.

30813. Cascales E, Lloubes R, Sturgis JN. (2001) The TolQ-TolR proteins energize TolA and share 309homologies with the flagellar motor proteins MotA-MotB. *Mol Microbiol* 42, 795-807.

31014. Celia H, Noinaj N, Zakharov SD, Bordignon E, Botos I, Santamaria M, Barnard TJ, Cramer 311WA, Lloubes R, Buchanan SK. (2016) Structural insight into the role of the Ton complex in 312energy transduction. *Nature* 538, 60-65.

31315. Pramanik A, Zhang F, Schwarz H, Schreiber F, Braun V. (2010) ExbB protein in the314cytoplasmic membrane of *Escherichia coli* forms a stable oligomer. *Biochemistry* 49, 8721-8728.

31516. Pramanik A, Hauf W, Hoffmann J, Cernescu M, Brutschy B, Braun V. (2011) Oligomeric
316structure of ExbB and ExbB-ExbD isolated from *Escherichia coli* as revealed by LILBID mass
317spectrometry. *Biochemistry* 50, 8950-8956.

31817. Sverzhinsky A, Fabre L, Cottreau AL, Biot-Pelletier DM, Khalil S, Bostina M, Rouiller I,
319Coulton JW. (2014) Coordinated rearrangements between cytoplasmic and periplasmic domains
3200f the membrane protein complex ExbB-ExbD of *Escherichia coli*. *Structure* 22, 791-797.

32118. Yonekura K, Maki-Yonekura S, Homma M. (2011) Structure of the flagellar motor protein 322complex PomAB: implications for the torque-generating conformation. *J Bacteriol* 193, 3863-3233870.

32419. Kim JS, Jeong H, Song S, Kim HY, Lee K, Hyun J, Ha NC. (2015) Structure of the tripartite 325multidrug efflux pump AcrAB-TolC suggests an alternative assembly mode. *Mol Cells* 38, 180-326186.

32720. Low HH, Gubellini F, Rivera-Calzada A, Braun N, Connery S, Dujeancourt A, Lu F, Redzej 328A, Fronzes R, Orlova EV, Waksman G. (2014) Structure of a type IV secretion system. *Nature* 329508, 550-553.

33021. Durand E, Nguyen VS, Zoued A, Logger L, Péhau-Arnaudet G, Aschtgen MS, Spinelli S,
331Desmyter A, Bardiaux B, Dujeancourt A, Roussel A, Cambillau C, Cascales E, Fronzes R.
332(2015) Biogenesis and structure of a type VI secretion membrane core complex. *Nature* 523,
333555-560.

33422. Bakelar J, Buchanan SK, Noinaj N. (2016) The structure of the β-barrel assembly machinery 335complex. *Science* 351, 180-186.

33623. Unger T, Jacobovitch Y, Dantes A, Bernheim R, Peleg Y. (2010) Applications of the
337Restriction Free (RF) cloning procedure for molecular manipulations and protein expression. *J*338*Struct Biol* 172, 34-44.

33924. Tabor S, Richardson CC. (1985) A bacteriophage T7 RNA polymerase/promoter system for 340controlled exclusive expression of specific genes. *Proc Natl Acad Sci U S A*. 82, 1074-1078.

34125. Studier FW, Moffatt BA. (1986) Use of bacteriophage T7 RNA polymerase to direct 342selective high-level expression of cloned genes. *J Mol Biol* 189, 113-130.

34326. Narayanan A, Ridilla M, Yernool DA. (2011) Restrained expression, a method to
3440verproduce toxic membrane proteins by exploiting operator-repressor interactions. *Protein Sci.*34520, 51-61.

34627. Stuchlík S, Turna J. (1998) Overexpression of the FNR protein of *Escherichia coli* with T7 347expression system. *Folia Microbiol (Praha)* 43, 601-604.

34828. Doherty AJ, Connolly BA, Worrall AF. (1993) Overproduction of the toxic protein, bovine
349pancreatic DNaseI, in *Escherichia coli* using a tightly controlled T7-promoter-based vector.
350*Gene* 136, 337-340.

35129. van den Ent F, Löwe J. (2006) RF cloning: a restriction-free method for inserting target 352genes into plasmids. *J Biochem Biophys Methods* 67, 67-74.

35330. Jeong J, Yim H, Ryu J, Lee H, Lee J, Seen D, Kang S. (2012) One-step sequence- and 354ligation-independent cloning as a rapid and versatile cloning method for functional genomics 355studies. *Appl Environ Microbiol* 78, 5440-5443.

35631. Bénédetti H, Lazdunski C, Lloubes R. (1991) Protein import into Escherichia coli: colicins A 357and E1 interact with a component of their translocation system. *EMBO J* 10, 1989-1995.

35832. Derouiche R, Bénédetti H, Lazzaroni JC, Lazdunski C, Lloubes R. (1995) Protein complex 359within *Escherichia coli* inner membrane. TolA N-terminal domain interacts with TolQ and TolR 360proteins. *J Biol Chem* 270, 11078-11084.

36133. Konovalova A, Silhavy TJ. (2015) Outer membrane lipoprotein biogenesis: Lol is not the 362end. *Philos Trans R Soc Lond B Biol Sci* 370, 1679.

36334. Ringquist S, Shinedling S, Barrick D, Green L, Binkley J, Stormo GD, Gold L. (1992)
364Translation initiation in *Escherichia coli:* sequences within the ribosome-binding site. *Mol*365*Microbiol* 6, 1219-1229.

366

#### 367i. Acknowledgements

368We would like to thank M. Petiti, L. Houot, H. Célia and D. Duché for careful reading. ED and 369RL are are funded by the Centre National de la Recherche Scientifique, the Aix-Marseille 370Université and two grants from the Agence Nationale de la Recherche (ANR-10-JCJC-1303-03 371and ANR-14-CE09-0023, respectively)



372

373Figure 1: Cloning of the synthetic *tssJLM* operon encoding the T6SS membrane core

374**complex.** (**a**) The *tssL*, *tssM* and *tssJ* genes were PCR amplified from the entero-aggregative *sci*-375*I* T6SS operon. DNA sequences were added encoding an optimized RBS at the 5' end of each 376gene and sequences encoding a C-terminal Strep-II tag for TssJ and N-terminal 6×His and Flag 377tags for TssM and TssL, respectively. (**b**) The *tssJLM* artificial operon was cloned in the pRSF-378Duet vector. (**c**) Schemes representing the protein constructs, the sub-domain boundaries and 379some indicated characteristics (TM, transmembrane segments; SP, signal peptide; CYS, acylated 380cysteine).

381