

Isolation and culture of Methanobrevibacter smithii by co-culture with hydrogen-producing bacteria on agar plates

S.I. Traore, S. Khelaifia, N. Armstrong, J.C. Lagier, Didier Raoult

► To cite this version:

S.I. Traore, S. Khelaifia, N. Armstrong, J.C. Lagier, Didier Raoult. Isolation and culture of Methanobrevibacter smithii by co-culture with hydrogen-producing bacteria on agar plates. Clinical Microbiology and Infection, 2019, 10.1016/j.cmi.2019.04.008. hal-02202491

HAL Id: hal-02202491 https://amu.hal.science/hal-02202491

Submitted on 21 Dec 2021

HAL is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers. L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.



Distributed under a Creative Commons Attribution - NonCommercial 4.0 International License

1	Isolation and culture of Methanobrevibacter smithii by co-culture with hydrogen-
2	producing bacteria on agar plates
3	Sory Ibrahima TRAORE ¹ , Saber KHELAIFIA ¹ , Nicholas ARMSTRONG ¹ , Jean Christophe
4	LAGIER ¹ , Didier RAOULT ^{1*}
5	
6	1. Aix Marseille Univ, IRD, APHM, MEPHI, IHU-Méditerranée Infection, 19-21
7	Boulevard Jean Moulin, 13385, Marseille cedex 05, France
8	
9	*Corresponding author: Prof. Didier RAOULT, Aix Marseille Univ, IRD, APHM, MEPHI,
10	IHU-Méditerranée Infection 19-21 Bd Jean Moulin 13005 Marseille, France. Tel: +33 (0)4 13
11	73 24 01 fax: + 33 (0) 13 73 24 02.
12	Email: didier.raoult@gmail.com
13	
14	Keywords: Methanobrevibacter smithii, anaerobes, co-culture, hydrogen production, methane
15	production.
16	
17	Word account: 1884
18	Abstract word account: 229
19	Figure: 03

20 Abstract

21 **Objective**

22 Methanogenic archaea are considered as extremely oxygen-sensitive organisms and their

23 culture is fastidious, requiring specific equipment. We report here the conditions allowing the

24 cultivation of *M. smithii* in an anaerobic chamber without hydrogen addition.

25 Methods

We first enriched the stool sample in an anaerobic liquid medium, and then, to cultivate *M. smithii* with *B. thetaiotaomicron* and other hydrogen producing bacteria as, on solid medium in an anaerobic chamber, we divided the agar plats in two compartments and we seeded each strain on each compartment. Methane production was assessed by gas chromatographic and the growing colonies were authenticated by MALDI-TOF MS.

31 **Results**

We successfully cultured *Methanobrevibacter smithii* from a liquid culture medium inoculated with a stool collected from a healthy donor under an anaerobic chamber. The isolation in pure culture was then a success on agar medium by performing a co-culture with *Bacteroides thetaiotaomicron*. We have also successfully tested the co-cultivation of *M*. *smithii* with other known hydrogen-producing bacteria. Gas chromatographic tests showed that these strains produced hydrogen in different amounts and agar colonies of methanogens were obtained by co-culture with these bacteria and methane production was detected.

39 Conclusions

In this work, we propose a new approach to isolate and cultivate new strains of *M. smithii* by
using a co-culture based technique that can facilitate and make available the isolation of new
methanogenic archaea strains in clinical microbiology laboratories.

43 Introduction

The cultivation of methanogens is considered difficult and has required until recently 44 the presence of highly specific equipment; both a source of hydrogen and the evacuation of 45 oxygen (1). The absence of oxygen and the presence of hydrogen are required for the culture 46 of methanogens. Recently, we have been able to show that anaerobes are aero-intolerant by 47 48 culturing them in the presence of antioxidants that allowed culturing most of the anaerobic bacteria (2). The cultivation of some of these anaerobic bacteria generates hydrogen, 49 especially for *Bacteroides*, as previously described (3). This particularity allowed us to 50 develop a culture system for methanogens using a double chamber comprising a medium 51 52 containing antioxidants, allowing the culture of anaerobes (4). An antibiotic cocktail was added to the agar medium for growing Methanobrevibacter to inhibit methanogen-associated 53 bacteria in multi-microbial specimens such as stool specimens (4). This work was the subject 54 55 of a publication that reported a series of new methanogen isolates (4). We performed in our laboratory the analysis of gut microbiota by culturomics using anaerobic host (5). We were 56 57 recently surprised to discover by chance that one of the colonies identified on agar by MALDI-TOF MS (Matrix-Assisted Laser Desorption-Ionization Time-Of-Flight Mass 58 Spectrometry) was Methanobrevibacter smithii, without having taken special conditions for 59 its cultivation. This surprising discovery led us to believe that in this anaerobic chamber, 60 which also contained cultures of anaerobic bacteria, the presence of hydrogen had been 61 sufficient to allow the cultivation of *M. smithii*. In this work, we report the conditions that 62 could allow the cultivation of *M. smithii* in an ordinary anaerobic chamber without the 63 64 introduction of hydrogen.

3

65 Materials and Methods

66 Ethics and sample collection

A stool specimen was collected from a healthy 26-year-old Malian female living in France for the last three months. Prior to collection, the donor had given his informed consent and the study was validated by the Ethics Committee of the IHU Méditerranée Infection under number 2016-011. The fecal specimen was collected in a recipient including a small anaerobic generator bag GENbag anaer (bioMérieux) immediately introduced to eliminate oxygen. Thirty minutes after the sampling, the stool sample was transported to the anaerobic chamber where we proceeded to its culture.

74 Culture procedure

One gram of stool specimen was suspended in 2 mL of phosphate buffered saline (PBS) (Life 75 Technologies, Carlsbad, CA, USA) and inoculated in an anaerobic blood culture bottle 76 77 (Becton Dickinson, Pont de Claix, France) supplemented with 5% of a mircofiltered (0.2µm) rumen. After 72 hours of incubation at 37 ° C inside the anaerobic chamber, we carried out 10 78 79 decreasing serial dilutions using PBS, then 50µL were seeded on 5% sheep blood-enriched Columbia agar (bioMérieux, Marcy l'Etoile, France) (COS). Agar plates were then incubated 80 at 37 ° C inside of the anaerobic chamber (Don Whitley Scientific, France). After 4 days of 81 82 incubation, the growing colonies were subcultured onto sterile agar plates for purification. The COS agar plate was divided into 8 equal parts and each isolated colony was subcultured 83 onto one part of the agar (Figure 1) and then incubated at 37 ° C for 96 hours. 84

85 Co-culture of Methanobrevibacter smithii and Bacteroides thetaiotaomicron

86 In order to co-cultivate *M. smithii* with *B. thetaiotaomicron* on solid medium, we divided the

- agar plats in two compartments, then we seeded *M. smithii* on one compartment and *B*.
- thetaiotaomicron on the other. The agar plates were then incubated for 72 hours at $37 \degree C$ in a
- sealed plastic bag with an anaerobiosis generator GENbag anaer (bioMérieux). After 72 hours

- 90 of incubation, the gas atmosphere was analyzed by gas chromatography (Figure 1) and the
- 91 growing colonies were identified by Maldi TOF-MS as previously described (5).
- 92 We also co-cultured *M. smithii* with *B. thetaiotaomicron* in liquid medium; an anaerobic
- 93 blood culture flask containing *M. smithii* and *B. thetaiotaomicron* was prepared and incubated
- 94 at 37 $^{\circ}$ C for 5 days. Then, we proceeded to the measurement of hydrogen and methane
- 95 contained in the blood culture flasks by mass spectrometry. The anaerobic blood culture flask
- 96 inoculated with *M. smithii* alone was introduced as negative control.

97 Co-culture of *M. smithii* with hydrogen-producing bacteria

- 98 The co-culture of *M. smithii* with hydrogen-producing bacteria was reported in the anaerobic
- 99 chamber on 5% sheep blood-enriched Columbia agar medium with the following bacterial
- 100 species: Bacteroides fragilis, Bacteroides vulgatus, Parabacteroides distasonis,
- 101 Marseillibacter massiliensis, Allisonella histaminiformansas and Desulfovibrio piger. The
- agar plates, after seeding the strains, were sealed in a plastic bag with anaerobic generator
- 103 GENbag anaer (bioMérieux) and incubated at 37 ° C for 72 hours.

104 MALDI-TOF MS identification

- 105 MALDI-TOF MS (Matrix-Assisted Laser Desorption-Ionization Time-Of-Flight Mass
- 106 Spectrometry) identification was adapted for isolated colonies as previously described (6).
- 107 Isolated colonies were deposited in duplicate on a MALDI-TOF target for further analysis.
- 108 Then, a Matrix solution enabled the ionization and desorption of the archaeal or bacterial
- 109 colonies were deposited on it. A score above 1.9 was considered as proper identification (7).

110 Hydrogen and methane measurements

- 111 Hydrogen and methane were measured using a Clarus 580 gas chromatography system
- 112 (Perkin Elmer, Villebon-sur-Yvette, France). 100 µL of gas was sampled from Blood culture
- bottle or airtight bags with a gastight syringe, then directly injected (split 10 mL/min) onto a
- 114 Shincarbon ST 80/100 micropacked column (2 m x 0.53 mm; Restek, Lisses, France). Injector

and oven were maintained at 110 and 70 ° C respectively. Gases were eluted using either 115 116 Helium or Argon depending on the compound detected by a Thermal Conductivity Detector (TCD) maintained at 100 ° C. Hydrogen was detected in the negative polarity at 40 mA using 117 Argon as elution/reference gas, whereas methane was detected in the positive polarity at 160 118 mA using Helium. Data recording and processing was performed using Totalchrom 6.3.2 119 (Perkin Elmer). A standard gas mixture containing 25 % of hydrogen, nitrogen, methane and 120 121 carbon dioxyde (Air Products, Allentown, USA) was used to estimate the amounts of gases. Results 122

123 Isolation procedure

Here, a co-culture based procedure was used to isolate and cultivate *M. smithii* from the stoolspecimen of a 26-year-old healthy Malian female.

The anaerobic Blood culture bottle inoculated by the stool specimen was methane-positive 126 after 72 hours of incubation at 37 ° C. Then, after serial dilutions, colonies appeared on the 127 agar medium after four days of incubation at 37 ° C in the anaerobic chamber. The growing 128 129 colonies subcultured onto sterile agar plates were strewn with numerous microcolonies that could not be distinguished from each other with the naked eye (Figure 1). The purification 130 step (Figure 1) allowed the isolation of one archaeal species and different bacterial species 131 132 including: Methanobrevibacter smithii, Allisonella histaminiformans, Marseillibacter massiliensis, Bacteroides fragilis, Bacteroides vulgatus, Parabacteroides distasonis and 133 Desulfovibrio piger. These bacteria all grew on the same plate as M. smithii come from the 134 same sample. 135

136 Co-culture of Methanobrevibacter smithii and Bacteroides thetaiotaomicron

137 Here, *B. thetaiotaomicron*, was used as hydrogen generator required for the growth of *M*.

smithii as previously described (4). After 72 hours of incubation, colonies grew on the two

agar plates compartments. Analysis of the gas atmosphere by gas chromatography revealed

140 the presence of methane, characteristic of methanogenic archaea growth (Figure 1).

141 Identification by MALDI-TOF MS then confirmed the authenticity of the two cultivated

142 microorganisms. *M. smithii* colonies were about 0.5 μm in diameter and ocher (Figure 1).

143 Methane was detected in the blood culture bottle when *M. smithii* was co-cultured with *B*.

144 *thetaiotaomicron*. The quantity of detected methane increased if microfiltred rumen was

added (data not shown). No methane was detected when *M. smithii* was cultivated alone in a

146 blood culture bottle.

147 Co-culture of *M. smithii* with hydrogen-producing bacteria

In order to determine which strain was able to induce the growth of *M. smithii*, each of these 148 149 strains (Bacteroides fragilis, Bacteroides vulgatus, Parabacteroides distasonis, Desulfovibrio piger Marseillibacter massiliensis and Allisonella histaminiformans) was subcultured on agar 150 plate with *M. smithii* using the same method as the co-culture with *B. thetaiotaomicron*. Then, 151 152 the methane and hydrogen production was measured by gas chromatography and the authenticity of the growing colonies was confirmed by MALDI-TOF MS. Growth was 153 observed after three days of incubation at 37 ° C when M. smithii was co-cultivated with B. 154 fragilis, B. vulgatus, P. distasonis and D. piger. Methane was detected by gas 155 chromatography, and represented a production of approximately 50 to 100 ppm (Figure 2). 156 157 No growth was observed when *M. smithii* was co-cultured with *M. massiliensis* and *A.* histaminiformans. 158

159 **Discussion**

160 Here, we successfully cultivated *M. smithii* in liquid enrichment and then in pure culture on

solid medium without the need for an external supply of hydrogen. This is due to the

162 effectiveness of a culture technique based on a co-culture carried out between a hydrogen-

163 consuming methanogenic strain (*M. smithii*) and a hydrogen-producing bacterium (*B.*

164 *thetaiotaomicron*). *M. smithii* is known to colonize 97.4% of the human population (8). After

an expected culture of *M. smithii* on agar medium without any addition of hydrogen, we 165 166 successfully cultivated this methanogenic archaea in a liquid medium by culture enrichment of stool specimen collected from a healthy donor, to which we added B. thetaiotaomicron as 167 hydrogen-generating bacterium. This enrichment was methane-positive after three days of 168 incubation, indicating the growth of methanogenic archaea (9). The isolation in pure culture 169 170 was then successfully obtained on agar medium still in co-culture with *B. thetaiotaomicron*. 171 Subsequently, we were able to subculture this isolate in co-culture with other known hydrogen-producing bacteria isolated at the same time and on the same agar plate as M. 172 smithii. We then tested the effect of each of these strains on the growth of M. smithii by 173 174 performing co-culture on solid medium. Growth of M. smithii was observed with each of these strains and gas chromatography measurements showed that these strains produced 175 hydrogen. This gas was detected in different amounts ranging from 0 to 75 ppm (Figure 3). 176 177 Viability and growth were confirmed by observing colony formation on solid medium and methane production was measured gas chromatography. 178

179 Culture-based techniques for the identification of methanogenic archaea mainly use the 180 culture technique in Hungate tubes developed by Hungate et al 1969 (2). This technique remains fastidious and not accessible to some clinical microbiology laboratories that do not 181 182 have the necessary equipment for its application despite evidence of presence of archaea in pathologies such as brain abscesses, paravertebral muscle abscess or periodontis as recently 183 demonstrated (10-12). Recently, we have developed in our laboratory a new culture technique 184 in a double chamber, where methanogenic archaea are co-cultured with Bacteroides 185 thetaiotaomicron which produces the hydrogen necessary for their growth (4). This technique 186 subsequently allowed us to isolate a large number of *M. smithii* strains from different samples 187 collected from different human mucosae (oral cavity, gut and vagina). In this work, and 188 thanks to the acquisition of a new anaerobic chamber in our laboratory, we propose a new 189

8

- approach to isolate and cultivate *M. smithii* isolates. The simplicity of this technique and the
- 191 few means necessary to put it into practice make the cultivation of methanogenic archaea
- 192 within the reach of all clinical microbiology laboratories that do not have specific equipment
- 193 enabling them to study methanogenic archaea.


194 Acknowledgments

- 195 This work has benefited from the French State support, managed by the 'Agence Nationale
- 196 pour la Recherche' including the "Programme d'Investissement d'avenir" under the reference
- 197 Méditerranée Infection 10-IAHU-03.
- 198 This work was supported by Région Provence Alpes Côte d'Azur and European funding
- 199 FEDER PRIMI.
- 200 **Conflict of Interest**
- 201 There is no conflict of interest
- 202

203		Reference List
204		
205 206	1	Dridi B, Raoult D, Drancourt M. 2011. Archaea as emerging organisms in complex human microbiomes. Anaerobe. 17(2):56-63.
207 208 209	2	Hungate RE. Chapter IV A Roll Tube Method for cultivation of strict anaerobes. 1969. In: Norris JR, Ribbons DW, editors. Methods in Microbiology. 3 ed. Academic Press, p. 117-32.
210 211	3	Samuel BS, Gordon JI. 2006. A humanized gnotobiotic mouse model of host-archaeal- bacterial mutualism. Proc Natl Acad Sci U S A. 103(26):10011-6.
212 213 214	4	Khelaifia S, Lagier JC, Nkamga VD, Guilhot E, Drancourt M, Raoult D. 2016. Aerobic culture of methanogenic archaea without an external source of hydrogen. Eur J Clin Microbiol Infect Dis. 35(6):985-91.
215 216 217	5	Lagier JC, Dubourg G, Million M, Cadoret F, Bilen M, Fenollar F, Levasseur A, Rolain JM, Fournier PE, Raoult D. Culturing the human microbiota and culturomics. Nat Rev Microbiol. 2018 May 1:540-550. doi: 10.1038/s41579-018-0041-0. Review.
218 219 220 221	6	Seng P, Abat C, Rolain JM, Colson P, Lagier JC, Gouriet F, Fournier PE, Drancourt M, La Scola B, Raoult D. 2013. Identification of rare pathogenic bacteria in a clinical microbiology laboratory: impact of Matrix-Assisted Laser Desorption Ionization Time of Flight Mass Spectrometry. J Clin Microbiol. 51(7):2182-94.
222 223 224 225	7	Elsawi Z, Togo AH, Beye M, Dubourg G, Andrieu C, Armsrtong N, Richez M, di Pinto F, Bittar F, Labas N, Fournier PE, Raoult D, Khelaifia S. 2017. <i>Hugonella</i> <i>massiliensis</i> gen. nov., sp. nov., genome sequence, and description of a new strictly anaerobic bacterium isolated from the human gut. MicrobiologyOpen. 6(4):e00458.
226 227 228	8	Dridi B, Henry M, El Khechine A, Raoult D, Drancourt M. 2009. High prevalence of <i>Methanobrevibacter smithii</i> and <i>Methanosphaera stadtmanae</i> detected in the human gut using an improved DNA detection protocol. PLOS ONE. 4(9):e7063.
229 230	9	Khelaifia S, Raoult D, Drancourt M. 2013. A versatile medium for cultivating methanogenic archaea. PLOS ONE. 8(4):e61563.
231 232 233 234	10	Drancourt M, Nkamga VD, Lakhe NA, Régis JM, Dufour H, Fournier PE, Bechah Y, Scheld WM, Raoult D. 2017. Evidence of archaeal methanogens in brain abscess. Clin Infect Dis. 65(1):1-5.
235 236 237 238 239	11	Nkamga VD, Lotte R, Roger PM, Drancourt M, Ruimy R. <i>Methanobrevibacter smithii</i> and <i>Bacteroides thetaiotaomicron</i> cultivated from a chronic paravertebral muscle abscess. Clin Microbiol Infect. 2016 Dec;22(12):1008-1009.
239 240 241 242 243	12	Huynh HTT, Pignoly M, Drancourt M, Aboudharam G. A new methanogen " <i>Methanobrevibacter massiliense</i> " isolated in a case of severe periodontitis. BMC Res Notes. 2017 Dec 1;10(1):657.

244 Figures legends


- Figure 1: Illustration of the different steps of the isolation procedure and the cultivation of *M*. *smithii* by co-culturing with *B. thetaiotaomicron*.
- **Figure 2:** Methane measurements by gas chromatography. Each plot represents the gas phase
- from the co-culture of *M. smithii* with: A, *P. distasonis*; B, *B. fragilis*; C, *B. vulgates*; D, *D.*
- 249 *piger*.
- 250 Figure 3: Hydrogen measurements by gas chromatography. Each plot represents the gas
- 251 phase from the co-culture of *M. smithii* with: A, *P. distasonis*; B, *B. fragilis*; C, *B. vulgates*;
- 252 D, *D. piger*.




and B. thetaiotaomicron

Subculuring *M. smithii* by co-culture with *B. thetaiotaomicron*

MALDI-TOF MS-based identification





ICD intensity (µV)