

SCA Medium: A New Culture Medium for the Isolation of All Candida auris Clades

Ahmad Ibrahim, Lucie Peyclit, Rim Abdallah, Saber Khelaifia, Amanda Chamieh, Jean-Marc Rolain, Fadi Bittar

► To cite this version:

Ahmad Ibrahim, Lucie Peyclit, Rim Abdallah, Saber Khelaifia, Amanda Chamieh, et al.. SCA Medium: A New Culture Medium for the Isolation of All Candida auris Clades. Journal of Fungi, 2021, 7 (6), pp.433. 10.3390/jof7060433 . hal-03333972

HAL Id: hal-03333972 https://amu.hal.science/hal-03333972

Submitted on 24 Sep 2021

HAL is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers. L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.



Distributed under a Creative Commons Attribution 4.0 International License





Brief Report SCA Medium: A New Culture Medium for the Isolation of All Candida auris Clades

Ahmad Ibrahim ^{1,2}, Lucie Peyclit ^{1,2}, Rim Abdallah ^{1,2}, Saber Khelaifia ^{1,2}, Amanda Chamieh ^{1,2}, Jean-Marc Rolain ^{1,2} and Fadi Bittar ^{1,2,*}

- ¹ Aix Marseille University, IRD, APHM, MEPHI, 13005 Marseille, France; ahmad.ibrahim@etu.univ-amu.fr (A.I.); lucie.peyclit@ap-hm.fr (L.P.); rim.abdallah@etu.univ-amu.fr (R.A.); saber.khelaifia@univ-amu.fr (S.K.); amanda.chamieh@etu.univ-amu.fr (A.C.); jean-marc.rolain@univ-amu.fr (J.-M.R.)
- ² IHU-Méditerranée Infection, 13005 Marseille, France
- * Correspondence: fadi.bittar@univ-amu.fr

Abstract: *Candida auris* is an emerging multidrug-resistant yeast causing nosocomial infections and associated with high mortality in immunocompromised patients. Rapid identification and characterisation are necessary for diagnosis and containing its spread. In this study, we present a selective culture medium for all *C. auris* clades. This medium is sensitive with a limit of detection ranging between 10^1 and 10^2 CFU/mL. The 100% specificity of SCA (specific *C. auris*) medium is confirmed on a set of 135 *Candida* strains, 50 bacterial species and 200 human stool samples. Thus, this medium specifically selects for *C. auris* isolation from clinical samples, allowing the latter to study its phenotypic profile.

Keywords: Candida; Candida auris; culture; emerging fungus; isolation; specific medium

1. Introduction

Candida auris is an emerging multidrug-resistant pathogen that was first isolated in 2009 [1] and is now known to have four geographical clades: South Asia: India; East Asia: Japan; Southern Africa: South Africa; and South America: Venezuela [2]. Recently, a potential fifth clade (Clade V) has been described in Iran [3]. *C. auris* is a biofilm-forming, halotolerant, thermo-resistant yeast [1,4] that can grow at temperatures ranging between 30 °C and 42 °C and can tolerate up to 10% salinity [5]. This promotes the ability of *C. auris* to colonize various medical equipment, plastic surfaces, and nosocomial environments [2–4], rendering *C. auris* an invasive, powerful pathogen. In addition, *C. auris* is highly resistant to different classes of antifungal agents such as azoles, amphotericin B and echinocandins [6] and is associated with significant mortality, especially in immuno-compromised patients with multiple comorbidities such as diabetes mellitus, renal failure, and cardiovascular disease [2,4]. Thus, its rapid identification and characterization are necessary to optimize clinical outcomes and to attempt to contain its nosocomial and worldwide spread.

The diagnosis of *C. auris* is challenging. The vast majority of commercial, readily available diagnostics, such as VITEK, API20C-AUX, Auxa-Color 2, BD Phoenixand MicroScan are misleading. *C. auris* is misidentified as other *Candida* species, namely *Candida famata*, *Candida haemulonii*, *Candida duobushaemulonii*, *Candida sake*, *Candida lustaniae* and *Candida guilliermondii* among others [7,8]. This often leads to delays in appropriate management and treatment. However, the use of Matrix-Assisted Laser Desorption/Ionization Time of Flight Mass Spectrometry (MALDI-TOF-MS) is promising [7,8], since the spectrum of *C. auris* is available in the reference database. In addition, the user can carefully select, "hand-pick", and isolate the corresponding identified colonies from other colonies growing on a plate. Moreover, the molecular identification and characterization of *C. auris* is well



Citation: Ibrahim, A.; Peyclit, L.; Abdallah, R.; Khelaifia, S.; Chamieh, A.; Rolain, J.-M.; Bittar, F. SCA Medium: A New Culture Medium for the Isolation of All *Candida auris* Clades. *J. Fungi* **2021**, *7*, 433. https:// doi.org/10.3390/jof7060433

Academic Editors: Jacques F. Meis and Anuradha Chowdhary

Received: 22 April 2021 Accepted: 28 May 2021 Published: 29 May 2021

Publisher's Note: MDPI stays neutral with regard to jurisdictional claims in published maps and institutional affiliations.



Copyright: © 2021 by the authors. Licensee MDPI, Basel, Switzerland. This article is an open access article distributed under the terms and conditions of the Creative Commons Attribution (CC BY) license (https:// creativecommons.org/licenses/by/ 4.0/). developed [9], and the use of some real-time PCR-targeting *C. auris* and sequencing of the ITS region is promising [9–12]. In 2017, a modified Sabouraud medium was suggested by Welsh et al. for a specific isolation of *C. auris* [5]. However, some *Candida* isolates were still misidentified [7,8]. Thus, our aim was to develop a specific medium only selective for all *C. auris* clades from clinical samples.

2. Materials and Methods

We collected 135 fungal strains (29 different species) [13–15] (Table 1) and cultured them on the modified Sabouraud broth suggested by Welsh et al. that contains 20 g/L Mannitol as a carbon source (in order to inhibit *C. glabrata* growth) [5] in both liquid and solid-phase media. For the liquid phase, turbidity of the broth was measured by spectrophotometry at 24, 48, and 72 h of incubation at 40 °C. This served to control the growth of each tested strain. For the solid phase, we added autoclaved bacterial agar (15 g/L) and adjusted for a pH = 7. Each colony was streaked directly on the solid agar or the liquid inoculum.

Table 1. Sti	rains and samp	les tested on t	the SCA medium.
--------------	----------------	-----------------	-----------------

Туре	Strain/Sample	Nb	Source	Origin	Identification	qPCR <i>C. auris</i> (Ibrahim et al., 2021)
	Bacillus cereus	1	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
	Corynebacterium amycolatum	1	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
	Corynebacterium jeikeium	1	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
	Corynebacterium propinquum	1	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
	Corynebacterium striatum	1	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
Gram-positive bacteria	Enterococcus faecalis	1	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
	Enterococcus faecium	1	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
	Micrococcus luteus	1	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
	Staphylococcus aureus	1	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
	Staphylococcus capitis Staphylococcus cohnii	1	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
		1	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
	Staphylococcus epidermidis	1	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
	Staphylococcus haemolyticus	1	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
	Staphylococcus lugdunensis	1	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
	Staphylococcus pasteuri	1	Clinical	Marseille, France	MALDI-TOF- MS	Not tested

Туре	Strain/Sample	Nb	Source	Origin	Identification	qPCR <i>C. auris</i> (Ibrahim et al., 2021)
	Staphylococcus saprophyticus	1	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
	Staphylococcus simulans	1	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
	Staphylococcus warneri	1	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
	Streptococcus agalactiae	1	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
	Streptococcus dysgalactiae	1	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
	Streptococcus equinus	1	Clinical	Marseille, France	MALDI-TOF- MS MALDI TOF	Not tested
	Streptococcus mitis	1	Clinical	France	MALDI-TOF- MS	Not tested
	Streptococcus pneumoniae Stanbulococcus	1	Clinical	France Marcoille	MALDI-TOF- MS MALDI TOF	Not tested
	hominis	1	Clinical	France Marsoille	MALDI-TOF- MS MALDI TOF	Not tested
Subtotal	salivarius	1	Clinical	France	MALDI-TOI- MS	Not tested
Subtotal	Achromobacter	2.5		Marseille,	MALDI-TOF-	NT + + - 1
	xylosoxidans	1	Clinical	France	MS	Not tested
	Acinetobacter baumannii	1	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
	Bacteroides fragilis	1	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
	Citrobacter braakii	1	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
	Citrobacter freundii	1	Clinical	France	MALDI-TOF- MS	Not tested
Carana arrastina	Citrobacter koseri Enterobacter	1	Clinical	France Marcoille	MALDI-TOF- MS MALDI TOF	Not tested
bacteria	aerogenes	1	Clinical	France	MALDI-TOF- MS	Not tested
	Enterobacter asburiae	1	Clinical	France	MALDI-TOF- MS	Not tested
	Enterobacter cloacae	1	Clinical	France	MALDI-TOF- MS	Not tested
	Enterobacter kobeii	1	Clinical	France	MALDI-TOF- MS	Not tested
	Escherichia coli	1	Clinical	France	MALDI-TOF- MS	Not tested
	influenzae	1	Clinical	France	MALDI-TOF- MS	Not tested
	parainfluenzae	1	Clinical	France	MALDI-TOF- MS	Not tested
	Hafnia alvei	1	Clinical	France	MALDI-TOF- MS	Not tested
	Klebsiella oxytoca	1	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
	Klebsiella pneumoniae	1	Clinical	Marseille, France	MALDI-TOF- MS	Not tested

Table 1. Cont.

Туре	Strain/Sample	Nb	Source	Origin	Identification	qPCR <i>C. auris</i> (Ibrahim et al., 2021)
	Moraxella catarrhalis	1	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
	Morganella morganii	1	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
	Pasteurella multocida	1	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
	Proteus mirabilis	1	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
	Proteus vulgaris	1	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
	Providencia stuartii	1	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
	Pseudomonas aeruginosa	1	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
	Raoultella ornithinolytica	1	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
	Stenotrophomonas maltophilis	1	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
Subtotal		25				
	Candida albicans	73	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
	Candida glabrata	8	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
	Candida krusei	4	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
	Candida parapsilosis	6	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
	Candida lusitaniae	3	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
	Candida tropicalis	6	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
	Candida zelanoides	1	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
	Candida lipolytica	1	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
Yeast	Candida inconspicua	1	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
	Candida intermedia	1	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
	Candida guilliermondii	3	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
	Candida bracarensis	1	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
	Candida utilis	1	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
	Candida bovina	1	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
	Candida dubliniensis	2	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
	Candida norvegensis	1	Clinical	Marseille, France	MALDI-TOF- MS	Not tested

Table 1. Cont.

Туре	Strain/Sample	Nb	Source	Origin	Identification	qPCR <i>C. auris</i> (Ibrahim et al., 2021)
	Candida kefyr	2	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
	Candida beverwijkiae	1	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
	Cryptococcus diffluens	1	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
	Cryptococcus uniguttulatus	1	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
	Cryptococcus neoformans	2	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
	Saccharomyces cerevisiae	2	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
	Rhodotorula mucilaginosa	1	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
	Yarrowia lipolitica	1	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
	Candida haemulonii	1	Clinical	Netherlands	MALDI-TOF- MS	Not tested
	Candida duobushaemulonii	1	Clinical	Netherlands	MALDI-TOF- MS	Not tested
	Trichosporon asahii	1	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
	Kodamaea ohmeri	1	Clinical	Marseille, France	MALDI-TOF- MS	Not tested
	<i>Candida auris</i> Clade I	2	Clinical	Netherlands	MALDI-TOF- MS	Positive
	<i>Candida auris</i> Clade II	1	Clinical	DSMZ collection	MALDI-TOF- MS	Positive
	<i>Candida auris</i> Clade III	2	Clinical	Netherlands	MALDI-TOF- MS	Positive
	<i>Candida auris</i> Clade IV	2	Clinical	Netherlands	MALDI-TOF- MS	Positive
Subtotal		135				
Human samples	Stool samples	200	Clinical	Marseille, France		Negative
Total		385				

Table 1. Cont.

We then cultured all *C. auris* (n = 7) and *Candida tropicalis* (n = 6) (see results) strains on the different culture media used in our diagnostic laboratories. The media used for culturing yeasts were: Sabouraud Dextrose Agar, CHROMagar *Candida* and Buffered Charcoal Yeast Extract. The media used for culturing bacteria were: Chocolate agar, Tryptic Soy agar, Columbia agar, Mannitol Salt agar, MacConkey agar (BioMérieux, Marcy-l'Etoile, France) and LBJMR (Lucie Bardet and Jean-Marc Rolain) medium [16].

Finally, we cultured this set of 13 strains (i.e., *C. tropicalis* and *C. auris*) on different combinations of the Welsh et al. broth and MacConkey (MCK) agar, a selective medium for Gram-negative bacteria (see results and discussion) (Table 2) [17]. We prepared several media with varying concentrations of bile salts (1.5 g/L, 1 g/L, 0.75 g/L) and/or crystal violet (0.5 mg/L) (used as inhibitors in MCK) [17] to the initial broth composition of Welsh et al. in solid and liquid phase with shaking (300 rpm) at 40 °C (Table 2).

Crystal Violet	Bile Salts	<i>C. auris</i> Clade I	<i>C. auris</i> Clade II	<i>C. auris</i> Clade III	<i>C. auris</i> Clade VI	C. tropicalis
0 mg/L	0 g/L	++	++	++	++	++
0 mg/L	0.5 g/L	++	++	++	++	++
0 mg/L	1 g/L	++	++	++	++	++
0 mg/L	1.5 g/L	++	+	++	++	+
0.5 mg/L	0 g/L	++	++	++	++	-
0.5 mg/L	0.5 g/L	++	++	++	++	-
0.5 mg/L	1 g/L	++	+	++	++	-
0.5 mg/L	1.5 g/L	++	+	++	++	-

Table 2. Growth results of *C. auris* (seven strains from four clades) and *C. tropicalis* (n = 6) tested strains according to crystal violet and bile salts concentrations for each condition 5 g of pancreatic digest of casein. 5 g of peptic digest of animal tissue, 100 g NaCl, 20 g of Mannitol, 50 mg/L chloramphenicol and 50 mg/L gentamicin were added.

To study the specificity of our designated medium, we cultured for 3 days at 40 °C a panel of various Gram-positive and Gram-negative bacteria (n = 50), different *C. auris* clades (DSM 21092 and 6 strains that were kindly provided by Dr. Jacques F Meis, Canisius Wilhemina Hospital, Department of Medical Microbiology and Infectious Diseases, the Netherlands) [9,18–20]. All *C. auris* strains were isolated from blood cultures except JCM 15,448 (Clade II), which is from the external ear canal. We tested also other *Candida* strains (including the closest *Candida* species to *C. auris* (*C. haemulonii* and *C. duobushaemulonii*) (directly streaked on the agar) [2,4]. In addition, we tested 200 fecal samples that were negative by real-time PCR for *C. auris*. These samples were collected at the Marseille Hospital (AP-HM, Assistance Publique-Hôpitaux Marseille) from routine laboratory diagnostics (Table 1). Given that this work did not involve human body/tissues or use clinical data from patients, and according to French law (Loi no 2012–300 of 5 March 2012 and Décret no 2016–1537 of 16 November 2016 published in the 'Journal Officiel de la République Française'), neither institutional ethical approval nor individual patient consent was required for this non-invasive study.

All patients routinely sign an approval that the samples they submit for testing at AP-HM labs may be used in research. Each stool sample was enriched in Tryptic soy broth (TSB) for 3 days at 37 °C. Then, a total of 0.1 mL of each enriched media was streaked on the solid media. All bacterial and fungal strains used in this work were selected from previous studies and identified correctly by MALDI-TOF-MS [13–15,21].

Moreover, to determine the sensitivity of our *C. auris*-specific medium, a series of ten-fold dilution $(10^{-1} \text{ to } 10^{-10})$ of 0.5 McFarland of each *C. auris* clade was cultured. For the second clade, we tested the strain of DSMZ collection: DSM 21092. In addition, we also cultured a mixture of the same ten-fold dilutions of 0.5 McFarland of each *C. auris* strain with a stool sample negative in real-time PCR for *C. auris* $(10^{-1} \text{ to } 10^{-10})$. We did this as an attempt to determine whether the natural presence of microbes in stools would affect the growth of *C. auris* on our designated medium.

3. Results

We first cultured all our available *Candida* spp. strains (n = 135) on the broth developed by Welsh et al. Interestingly, we isolated 100% of the tested *C. auris* (n = 7) and *C. tropicalis* (n = 6) species on the Welsh broth [5].

However, when we cultured all strains on a combination of MacConkey agar and the Welsh broth, the growth of *C. tropicalis* strains was totally inhibited. On the other hand, we observed a normal growth of all *C. auris* species.

The addition of crystal violet, with or without bile salts, did not affect the growth of *C. auris* (Figure 1), which we evaluated by measuring the turbidity of the broth using a spectrophotometer. However, we observed a total inhibition of all *C. tropicalis* strains in the presence of 0.5 mg/L crystal violet.

7 of 10



Figure 1. Striking *Candida auris* strains onto SCA medium, 2 strains/clade except for *C. auris* clade II (one strain).

The addition of bile salts only partially inhibited the growth of *C. auris* clade II. This was reflected as a decrease in the spectrophotometer value in comparison to the negative control only containing broth and crystal violet. As for *C. tropicalis* strains, the presence of bile salts without crystal violet did not inhibit their growth. It was then excluded from the final composition of the SCA medium (Table 2).

Therefore, the growth of *C. tropicalis* was inhibited by adding crystal violet at 0.5 mg/L to the initial broth suggested by Welsh et al. with no effect on the growth of the 7 *C. auris* strains (2 strains/clade: except for clade II (one strain)) in solid and liquid phase.

Moreover, the specificity of this medium was also confirmed after cultivating all bacterial spp. *Candida* spp. and fecal samples mentioned in Table 1 on the last composition of our medium (Table 3). No growth was noted after 3 days of incubation at 40 °C.

Table 3. Final composition of SCA (Specific Candida auris) medium.

	Pancreatic Digest of Casein	Peptic Digest of Animal Tissue	NaCl	Mannitol	Crystal Violet	Agar	рН	Chloramphenicol	Gentamicin
Welsh et al. broth	5 g	5 g	100 g	20 g	-	-	5.6	50 mg/L	50 mg/L
SCA medium	5 g	5 g	100 g	20 g	0.5 mg	15 g	7	50 mg/L	50 mg/L

Concerning the sensitivity of our medium, after cultivation of several dilutions of a solution of 0.5 McFarland (10^{-1} to 10^{-10}), the limit of detection (LOD) of *C. auris* Clade I, II and VI was 10^2 CFU/mL in both serial dilutions (with physiological water and a fecal sample). For the third clade (South Africa), we observed an excessive growth of both strains at an LOD of 10^1 CFU/mL. Thus, the presence of other microbes in a fecal sample does not inhibit or affect the growth of *C. auris*. All resulting colonies were isolated and subsequently identified by MALDI-TOF-MS and real-time PCR [9].

The final composition of our medium in 1 L of deionized water was 5 g of pancreatic digest of casein, 5 g of peptic digest of animal tissue, 20 g of mannitol, 0.5 mg of crystal violet (Sigma-aldrich, Darmstadt, Germany), 100 g NaCl, 50 mg/L of chloramphenicol and 50 mg/L of gentamicin with pH = 7 (\pm 0.2) at 40 °C (Figure 1, Table 3).

4. Discussion

Correct identification and rapid isolation of *C. auris* is essential for the timely and appropriate antifungal treatment and for infection prevention and control measures. This ensures limiting and controlling a possible nosocomial outbreak of a pathogen with fatal consequences and limited therapeutic options. This need is emphasized in the current COVID-19 pandemic we are living in, whereby an increasing amount of critically ill, immunocompromised patients populate the hospitals [22–24]. The gold standard of diagnosis remains molecular identification by sequencing and/or real-time PCR [9–12]. However, this is expensive and may only be available at specialized referral or research laboratories. Therefore, we have developed a cheap, simple, easily-prepared medium for isolating *C. auris* that is specific and sensitive.

A caveat to the enrichment broth developed by Welsh et al. for *C. auris* isolation is that it also grows *C. tropicalis* strains, which we demonstrated in our work. We show that the addition of crystal violet (0.5 mg/L) to our specific *C. auris* medium inhibited the growth of *C. tropicalis*. SCA (specific *C. auris*) is considered as a new version of the above-mentioned broth [5].

In addition, other selective media are recently described for the isolation of *C. auris* [25–27]. Generally, the selected criteria to cultivate this yeast remain the same as those used in our work: the thermo-resistance, halo-tolerance and multi-resistance to many antifungal agents [2,4]. Here, we managed to develop a medium with the minimum possible inhibitors (40 °C incubation instead of 42, and 10% salinity instead of 12.5% compared to another study [26]). These moderate inhibitors still maintain the selectivity of SCA medium for all *C. auris*.

Usually, *C. auris*, *C. krusei* and *C. parapsilosis* appear as pink colonies on the CHROMagar medium [2,4,28], and a newer version of CHROMagar has been developed to identify *C. auris* [25,27]. This phenotypic identification may be biased and reader-dependent, since it is based on interpretation of the color/morphology of growing colonies. Moreover, a longer duration of incubation may lead to a change in results as well. Thus, a high percentage of error may occur. Our SCA medium facilitates the interpretation of results by only selecting for *C. auris*.

Interestingly, spiking a stool sample with a *C. auris* strain affected neither the specificity nor the sensitivity of this medium. However, increasing our sample size and testing a larger amount of clinical *C. auris* isolates and related yeasts (such as *C. famata, C. pseudohaemulonii, C. metapsilosis, C. orthopsilosis, C. rugosa, C. vulturna, Wickerhamomyces anomalus* and *C. sake*) is necessary for further evaluation, validation and reproducibility of our *C. auris*specific medium.

5. Conclusions

The use of our developed medium enables the rapid, specific isolation of *C. auris* strains and helps in the timely management of patients and resources to limit the occurrence of *C. auris* outbreaks. We propose an implementation of SCA medium in routine clinical mycology for screening skin, urine, vaginal and blood samples, especially in high-risk populations [13,15]. This SCA medium will further enhance our understanding of the phenotypic characteristics of *C. auris* and future isolates, most importantly allowing an accurate study of their antifungal resistance profiles.

Author Contributions: Conceptualization, A.I., L.P., J.-M.R. and F.B.; methodology, A.I., R.A., S.K. and A.C.; writing—original draft preparation, A.I., L.P. and F.B.; writing—review and editing, J.-M.R. and F.B.; supervision, J.-M.R. and F.B. All authors have read and agreed to the published version of the manuscript.

Funding: This work was supported by the French Government under the Investissements d'avenir (Investments for the Future) program managed by the Agence Nationale de la Recherche (ANR, fr: National Agency for Research), (reference: Méditerranée Infection 10-IAHU-03). This work was supported by Région Provence-Alpes-Côte d'Azur and European funding (FEDER (Fonds européen de développement régional) PRIMMI (Plateformes de Rechercheet d'Innovation Mutualisées Méditerranée Infection)).

Institutional Review Board Statement: Ethical review and approval were waived for this study, as this research does not involve human body/tissue, and only discarded stool samples, once laboratory analyses had been initiated, were used. According to French law (Loi no 2012–300 du 5 March 2012 and Décret no 2016–1537 du 16 November 2016 published in the 'Journal Officiel de la République Française'), no Institutional Review Board is necessary in this case.

Informed Consent Statement: Patient consent is not required for this type of study (Loi no 2012–300 of 5 March 2012 and Décret no 2016–1537 of 16 November 2016 published in the 'Journal Officiel de la République Française').

Data Availability Statement: Not applicable.

Acknowledgments: We want to thank CookieTrad for English correction. We would like to thank Jacques F. Meis for sending us the strains of *C. auris*, *C. haemulonii* and *C. duobushaemulonii*.

Conflicts of Interest: The authors declare no conflict of interest.

References

- 1. Satoh, K.; Makimura, K.; Hasumi, Y.; Nishiyama, Y.; Uchida, K.; Yamaguchi, H. *Candida auris* sp. nov., a novel ascomycetous yeast isolated from the external ear canal of an inpatient in a Japanese hospital. *Microbiol. Immunol.* **2009**, *53*, 41–44. [CrossRef]
- Osei Sekyere, J. Candida auris: A systematic review and meta-analysis of current updates on an emerging multidrug-resistant pathogen. Microbiologyopen 2018, 7, e00578. [CrossRef]
- Chow, N.A.; De Groot, T.; Badali, H.; Abastabar, M.; Chiller, T.M.; Meis, J.F. Potential fifth clade of *Candida auris*, Iran, 2018. *Emerg. Infect. Dis.* 2019, 25, 1780–1781. [CrossRef]
- Jeffery-Smith, A.; Taori, S.K.; Schelenz, S.; Jeffery, K.; Johnson, E.M.; Borman, A.; Manuel, R.; Browna, C.S. Candida auris: A review of the literature. Clin. Microbiol. Rev. 2018, 31, e00029-17. [CrossRef] [PubMed]
- Welsh, R.M.; Bentz, M.L.; Shams, A.; Houston, H.; Lyons, A.; Rose, L.J.; Litvintseva, A.P. Survival, persistence, and isolation of the emerging multidrug-resistant pathogenic yeast *Candida auris* on a plastic health care surface. *J. Clin. Microbiol.* 2017, 55, 2996–3005. [CrossRef] [PubMed]
- Rudramurthy, S.M.; Chakrabarti, A.; Paul, R.A.; Sood, P.; Kaur, H.; Capoor, M.R.; Kindo, A.J.; Marak, R.S.K.; Arora, A.; Sardana, R.; et al. *Candida auris* candidaemia in Indian ICUs: Analysis of risk factors. *J. Antimicrob. Chemother.* 2017, 72, 1794–1801. [CrossRef] [PubMed]
- Kathuria, S.; Singh, P.K.; Sharma, C.; Prakash, A.; Masih, A.; Kumar, A.; Meis, J.F.; Chowdhary, A. Multidrug-resistant *Candida* auris misidentified as *Candida haemulonii*: Characterization by matrix-assisted laser desorption ionization-time of flight mass spectrometry and DNA sequencing and its antifungal susceptibility profile variability by vitek 2, CLSI broth microdilution, and etest method. J. Clin. Microbiol. 2015, 53, 1823–1830. [CrossRef] [PubMed]
- Mizusawa, M.; Miller, H.; Green, R.; Lee, R.; Durante, M.; Perkins, R.; Hewitt, C.; Simner, P.J.; Carroll, K.C.; Hayden, R.T.; et al. Can multidrug-resistant *candida auris* be reliably identified in clinical microbiology laboratories? *J. Clin. Microbiol.* 2017, 55, 638–640. [CrossRef] [PubMed]
- 9. Ibrahim, A.; Baron, S.A.; Yousfi, H.; Hadjadj, L.; Lalaoui, R.; Morand, S.; Rolain, J.M.; Bittar, F. Development and standardization of a specific real-time PCR assay for the rapid detection of *Candida auris*. *Eur. J. Clin. Microbiol. Infect. Dis.* **2021**, 1–5. [CrossRef]
- 10. Lima, A.; Widen, R.; Vestal, G.; Uy, D.; Silbert, S. A TaqMan probe-based real-time PCR assay for the rapid identification of the emerging multidrug-resistant pathogen candida auris on the BD max system. *J. Clin. Microbiol.* **2019**, 57. [CrossRef]
- Leach, L.; Zhu, Y.; Chaturvedi, S. Development and Validation of a Real-Time PCR Assay for Rapid Detection of *Candida auris* from Surveillance Samples. J. Clin. Microbiol. 2018, 56. [CrossRef] [PubMed]
- 12. Ahmad, A.; Spencer, J.E.; Lockhart, S.R.; Singleton, S.; Petway, D.J.; Bagarozzi, D.A.; Herzegh, O.T. A high-throughput and rapid method for accurate identification of emerging multidrug-resistant *Candida auris*. *Mycoses* **2019**, *62*, 513–518. [CrossRef] [PubMed]
- 13. Yousfi, H.; Ranque, S.; Rolain, J.M.; Bittar, F. In vitro polymyxin activity against clinical multidrug-resistant fungi. *Antimicrob. Resist. Infect. Control* **2019**, *8*. [CrossRef] [PubMed]
- 14. Yousfi, H.; Cassagne, C.; Ranque, S.; Rolain, J.M.; Bittar, F. Repurposing of ribavirin as an adjunct therapy against invasive Candida strains in an in vitro study. *Antimicrob. Agents Chemother.* **2019**, *63*. [CrossRef]
- 15. Lagier, J.C.; Khelaifia, S.; Alou, M.T.; Ndongo, S.; Dione, N.; Hugon, P.; Caputo, A.; Cadoret, F.; Traore, S.I.; Seck, E.H.; et al. Culture of previously uncultured members of the human gut microbiota by culturomics. *Nat. Microbiol.* **2016**, *1*, 16203. [CrossRef]
- 16. Bardet, L.; Le Page, S.; Leangapichart, T.; Rolain, J.M. LBJMR medium: A new polyvalent culture medium for isolating and selecting vancomycin and colistin-resistant bacteria. *BMC Microbiol.* **2017**, *17*, 220. [CrossRef]
- 17. Jung, B.; Hoilat, G.J. MacConkey Medium; n: StatPearls [Internet]; StatPearls Publishing: Treasure Island, FL, USA, 2021.
- 18. Calvo, B.; Melo, A.S.A.; Perozo-Mena, A.; Hernandez, M.; Francisco, E.C.; Hagen, F.; Meis, J.F.; Colombo, A.L. First report of *Candida auris* in America: Clinical and microbiological aspects of 18 episodes of candidemia. *J. Infect.* **2016**, *73*, 369–374. [CrossRef]
- Ruiz-Gaitán, A.; Moret, A.M.; Tasias-Pitarch, M.; Aleixandre-López, A.I.; Martínez-Morel, H.; Calabuig, E.; Salavert-Lletí, M.; Ramírez, P.; López-Hontangas, J.L.; Hagen, F.; et al. An outbreak due to *Candida auris* with prolonged colonisation and candidaemia in a tertiary care European hospital. *Mycoses* 2018, *61*, 498–505. [CrossRef]
- 20. Chowdhary, A.; Sharma, C.; Duggal, S.; Agarwal, K.; Prakash, A.; Singh, P.K.; Jain, S.; Kathuria, S.; Randhawa, H.S.; Hagen, F.; et al. New clonal strain of *Candida auris*, Delhi, India. *Emerg. Infect. Dis.* **2013**, *19*, 1670–1673. [CrossRef]
- L'Ollivier, C.; Cassagne, C.; Normand, A.C.; Bouchara, J.P.; Contet-Audonneau, N.; Hendrickx, M.; Fourquet, P.; Coulibaly, O.; Piarroux, R.; Ranque, S. A MALDI-TOF MS procedure for clinical dermatophyte species identification in the routine laboratory. *Med. Mycol.* 2013, 51, 713–720. [CrossRef]
- De Almeida, J.N.; Francisco, E.C.; Hagen, F.; Brandão, I.B.; Pereira, F.M.; Dias, P.H.P.; De Miranda Costa, M.M.; De Souza Jordão, R.T.; De Groot, T.; Colombo, A.L.; et al. Emergence of *Candida auris* in Brazil in a COVID-19 Intensive Care Unit. *J. Fungi* 2021, 7, 220. [CrossRef]
- 23. Allaw, F.; Kara Zahreddine, N.; Ibrahim, A.; Tannous, J.; Taleb, H.; Bizri, A.R.; Dbaibo, G.; Kanj, S.S. First *Candida auris* Outbreak during a COVID-19 Pandemic in a Tertiary-Care Center in Lebanon. *Pathogens* **2021**, *10*, 157. [CrossRef]

- Chowdhary, A.; Sharma, C.; Meis, J.F. Candida auris: A rapidly emerging cause of hospital-acquired multidrug-resistant fungal infections globally. PLoS Pathog. 2017, 13, e1006290. [CrossRef] [PubMed]
- De Jong, A.W.; Dieleman, C.; Carbia, M.; Tap, R.M.; Hagen, F. Performance of two novel chromogenic media for the identification of multidrug-resistant *candida auris* compared with other commercially available formulations. *J. Clin. Microbiol.* 2021, 59. [CrossRef]
- 26. Das, S.; Singh, S.; Tawde, Y.; Chakrabarti, A.; Rudramurthy, S.M.; Kaur, H.; Shankarnarayan, S.A.; Ghosh, A. A selective medium for isolation and detection of *Candida auris*, an emerging pathogen. *J. Clin. Microbiol.* **2021**, *59*. [CrossRef] [PubMed]
- 27. Borman, A.M.; Fraser, M.; Johnson, E.M. CHROMagarTM *Candida* Plus: A novel chromogenic agar that permits the rapid identification of *Candida auris*. *Med. Mycol.* **2021**, *59*, 253–258. [CrossRef] [PubMed]
- 28. Nadeem, S.G.; Hakim, S.T.; Kazmi, S.U. Use of CHROMagar *Candida* for the presumptive identification of *Candida* species directly from clinical specimens in resource-limited settings. *Libyan J. Med.* **2010**, *5*, 1–6. [CrossRef]